

— TRAPEZIOMETACARPAL JOINT PROSTHESIS

MAÏA™ *dual mobility*



DOCUMENTATION



groupe lépine

Depuis 1714

DUAL MOBILITY CUP

- titanium aluminum vanadium alloy.
- cementless fixation made possible by bilayer coating of vacuum-plasma sprayed pure porous titanium and hydroxyapatite.
- outer macrostructure: 4 polar stabilising spikes.
- sterilized by Gamma irradiation.



MAÏA™ DUAL MOBILITY CUP

Designation	Ref
MAIA DUAL MOBILITY CUP Ø9	M100 3009

DUAL MOBILITY LINER + NECK

- liner:
 - machined from UHMW-PE
- neck:
 - stainless steel
 - modular
 - offset respecting the anatomy
 - spherical 4-mm diameter head
- the liner and the neck are delivered preassembled.
- sterilized by Gamma irradiation.



MAÏA™ DUAL MOBILITY LINER + NECK

Designation	Ref
MAIA DUAL MOBILITY LINER M	M100 3190
MAIA DUAL MOBILITY LINER L	M100 3191
MAIA DUAL MOBILITY LINER XL	M100 3192

METACARPAL STEM

- titanium aluminum vanadium alloy.
- cementless fixation made possible by bilayer coating of vacuum-plasma sprayed pure porous titanium and hydroxyapatite.
- outer macrostructure: scales on longitudinal ridges providing stability.
- anatomically shaped design and homothetic range.
- sterilized by Gamma irradiation.



MAÏA™ METACARPAL STEM

Designation	Ref
MAIA STEM S7s	M100 0017
MAIA STEM S8s	M100 0018
MAIA STEM S9s	M100 0019
MAIA METACARPAL STEM S10	M100 0010

lépine Algérie
9 rue philosophe Tabrizi
Les sources
Bir Mourad Rais - ALGER
lepine-algerie@groupe-lepine.com

lépine Italia
Via Cassanese 100,
20090 Segrate (Milano)
lepine-italia@groupe-lepine.com

lépine Deutschland
Ohle Ring 23-25
Stade-Wiepenkathen, 21684
lepine-deutschland@groupe-lepine.com

lépine Maroc
79 avenue IBN SINA
10080 RABAT - AGDAL
lepine-maroc@groupe-lepine.com

lépine Ibérica
J.J. Tadeo Murguía, 3 y 5 Bajos
20304 IRÚN (GUIPÚZCOA)
lepine-iberica@groupe-lepine.com

175 rue Jacquard - CS 50307
69727 Genay Cedex - FRANCE
Tél. +33 (0)4 72 33 02 95
Fax +33 (0)4 72 35 96 50

www.groupe-lepine.com



groupe lépine

Depuis 1714